

**Manufacturer**

HEMODIA  
85 rue du chêne vert  
31670 Labège  
FRANCE

For more information :

France : +33 5 61 00 71 81

United Kingdom : +44 20 33 765 228

Germany : +49 69 80 92 94 71

Sweden : +46 313 04 50 23

Denmark : +45 787 877 05

Belgium : +32 28 08 56 11

Netherlands : +31 207 01 80 86


Switzerland : +41 435 50 99 09

HEMODIA is expanding and searching for more partnerships for Europe and International.

You wish to partner with HEMODIA ?

Please contact us on [hemodia@hemodia.com](mailto:hemodia@hemodia.com) or +33 561 00 71 81.

*Updated 11/16 - FT\_ART\_ANG\_2*



ARTHROSCOPIC  
TUBING



Technical brochure



## HEMODIA

Arthroscopy specialist since 1991

The arthroscopic tubing is manufactured and distributed by HEMODIA. HEMODIA was originally manufacturing and providing the Tubing used on FMS DUO® Pumps.


Hemodia ensures the maintenance of the FMS DUO®+ pump, manufactured by DePuySynthes, in accordance with the NF EN ISO 13485 quality standards..


DISTRIBUTEUR

FABRICANT


- 2014** HEMODIA integrates OXO MEDICAL and becomes independent by manufacturing and distributing his Tubing Portfolio.
- 2013** OXO MEDICAL launches his own tubing portfolio, manufactured by HEMODIA. OXO MEDICAL distributes the HEMODIA tubing all over Europe
- 2006** DepuyMitek acquires FMS and the HEMODIA Tubing contract
- 1996** Launch of DUO® + Pumps; tubing created and manufactured by HEMODIA
- 1991** Launch of FMS4® Pumps; tubing created and manufactured by HEMODIA









**DAY SET \*\***  
REF 45034504\_A (DEHP)  
REF 45034504\_B (DEHP free)  
30/box



**PATIENT SET WITH CHECK VALVE \***  
REF 45094549\_A (DEHP)  
REF 45094549\_B (DEHP free)  
30/box



**SINGLE SUCTION TUBE\***  
REF 45054505\_A (DEHP)  
REF 45054505\_B (DEHP free)  
40/box



**INTERMEDIARY TUBE WITH CHECK VALVE\*\***  
REF 11021142 \_A (DEHP)  
REF 11021142 \_B (DEHP free)  
50/box

HEMODIA Tubing is used with Pumps manufactured by DePuySynthes :  
 \*FMS DUO®+  
 \*\*FMS DUO®+ and FMS SOLO®

